


Jeffrey Meeusen¹, Frank Stämmler², John Lieske³, Marcello Grassi², Maulik Shah⁴, Laurence Dubourg⁵, and Eric Schiffer²

¹Mayo Clinic Department of Laboratory Medicine and Pathology (Rochester, MN), ²Department of Research and Development, numares AG (Regensburg, Germany), ³Mayo Clinic Division of Nephrology and Hypertension (Rochester, MN), ⁴Clinical Research, numares GROUP Corp. (Boston, MA), ⁵Service d'Explorations Fonctionnelles Rénales et Métaboliques, Hôpital Edouard Herriot (Lyon, France)

BACKGROUND

Accurate assessment of glomerular filtration rate (GFR) is critical to decision making in individuals with liver disease. Renal impairment is common in association with liver disease and the degree of renal dysfunction impacts decisions on drug dosing, therapeutic interventions, and suitability for liver transplantation. Altered hemodynamics in liver disease often results in overestimation of GFR when using creatinine based GFR estimating equations. Recently, we have developed a novel GFR equation (GFR_{NMR}), which utilizes serum myo-inositol, valine and creatinine quantified by nuclear magnetic resonance spectroscopy (NMR) in combination with Cystatin-C, age and sex. This equation outperforms many GFR estimating equations in chronic kidney disease.




METHODS

We compared various GFR estimation equations including GFR_{NMR} in end-stage liver disease patients scheduled for orthotopic liver transplantation (OLT) in a multicenter retrospective study. In, n=165, liver recipients, renal tracer clearance (mGFR by iothalamate, iohexol or inulin) was measured in preparation of OLT as part of normal clinical care to assess kidney function. We measured Cystatin-C and NMR kidney biomarkers, myo-inositol, valine and creatinine in a single simultaneous measurement. Analytes were used to estimate GFR based on CKD-EPI₂₀₀₉, CKD-EPI_{Cys}, or GFR_{NMR} .


	N
Age (yrs, range)	18 - 82
Age (yrs, mean ± SD)	58 ± 12
Sex (% male)	65.5
BMI range	16 - 52
BMI mean ± SD	29 ± 7
BMI: < 20	9 (5.5%)
BMI: 20-25	90 (54.5%)
BMI: > 30	66 (40%)
mGFR range	3 - 146
mGFR mean ± SD	65 ± 29
iothalamate	113 (68.5%)
iohexol	5 (3%)
inulin	47 (28.5%)
CKD Stage 1	34 (20.6%)
CKD Stage 2	52 (31.5%)
CKD Stage 3	60 (36.4%)
CKD Stage 4	18 (10.9%)
CKD Stage 5	1 (0.6%)

Characteristics of 165 patients with end-stage liver disease

RESULTS



GFR_{NMR} was superior to CKD-EPI₂₀₀₉ and CKD-EPI_{Cys} as shown by the percentage of samples with relative error < 15, 20, or 30% (P15, P20, P30).



GFR_{NMR} showed lower median bias (A) and higher accuracy (B, shown by 1-P20) than CKD-EPI₂₀₀₉ and CKD-EPI_{Cys} over the entire GFR range, especially in the clinically relevant range below 60 ml/min/1.73m².

Equation	range	mean±SD	P15 [95% CI]	P20 [95% CI]	P30 [95% CI]	MAE [95% CI]	median Bias [95% CI]
CKD-EPI ₂₀₀₉	8 - 142	72 ± 30	44.24 [36.36 - 52.12] ns	50.91 [43.03 - 58.79]*	65.45 [58.2 - 72.73]**	15.28 [12.94 - 17.35]*	3 [-1 - 6]***
CKD-EPI _{Cys}	9 - 120	52 ± 27	27.88 [21.21 - 34.55] ***	40 [33.33 - 47.27]***	61.21 [53.95 - 68.48]**	16.78 [14.83 - 18.58]***	-13 [-16 - -10]***
GFR_{NMR}	9 - 130	62 ± 26	53.33 [46.67 - 61.21]	64.85 [58.18 - 72.12]	80.61 [75.15 - 86.67]	11.23 [9.57 - 12.71]	-2 [-5 - 1]

Key Performance Indicators of GFR_{NMR} in comparison to CKD-EPI₂₀₀₉ and CKD-EPI_{Cys} (MAE: Mean absolute error; CI: Confidence interval)

Significance level adj. p-value, comparison with GFR_{NMR} : ns = not significant; * < 0.01; ** < 0.001; *** < 0.0001; significance of differences was calculated via McNemar test for P15, P20 and P30 or Wilcoxon-signed rank test for MAE and median Bias;

CONCLUSIONS

Kidney function assessment is critical in patients with Liver Diseases. GFR_{NMR} provides a more accurate assessment of GFR estimation over creatinine and Cystatin-C based equations in individuals with liver disease.

CONTACT DETAILS

maulik.shah@numares.com
www.numares.com

Follow numares: @numares AG
 @Agnumares

@numares AG
 numares NMR diagnostics

